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- 產品型號：1、BOD6060W □
2、BOD8060W □

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
User Manual ▶

Product Introduction

This product utilizes fiber optic micro-vibration sensing technology to achieve zero-contact and continuous monitoring of human vital signs data, providing continuous health monitoring, early disease prevention, and real-time alert services for people with chronic diseases such as cardiovascular and cerebrovascular diseases, snoring, post-surgery individuals, or elderly people living alone.

The product allows multiple monitoring pads to be managed through a single administrator account.

Main functions: respiratory monitoring, heart rate monitoring, in-bed/out-of-bed monitoring, and data anomaly alerts.

There are 4 indicator lights in the upper right corner of the product (from top to bottom):

- 1) **Power:** Solid green light indicates normal power supply.
- 2) **Error:** Solid red light indicates device malfunction.
- 3) **Connection:** Flashing green light indicates successful network pairing and normal connection to the server.
- 4) **Bluetooth:** Solid blue light indicates Bluetooth is enabled and network configuration is in progress.

Items included in the packaging box

Item	Quantity
Vital signs monitor	1 pcs
User Manual	1 pcs
Certificate of Conformity	1 pcs
Power Adapter	1 pcs

Product Specifications

Product Name	vital signs monitor
Model	BOD6060W
Size	L600*W600*H20mm
Material	Waterproof fabric, PVC frame, Glass fiber optic cable
Weight	2.5 kg
Maximum load capacity	150kg
Installation position	Placed under a mattress measuring 8cm-45cm
Suitable bed types	household beds or medical care beds
Storage temperature	-20°C - 85°C
Heart rate	Range: 40bpm~150bpm Error: ± 3 bpm or $\pm 3\%$,
Respiratory rate	Measurement range: 0-45rpm Measurement error: 6-45rpm: ± 2 rpm; 0-5rpm: Not pre-defined
In/out of bed monitoring	support
Body motion monitoring	support
Power supply	5V/2A power plug (adapter)
Network communication	Bluetooth 4.0 or above +WiFi (802.11/b/g/n-2.4G)

Product Usage Conditions

APP Download

Search for Anvenus Plus (Anvenus+) in the App Store or Google Play, then download and install it.

Installation method


Place monitoring mattresses

Remove the product (monitoring mattress) and place it between the mattress and the bed frame. Ensure the surface of the monitoring mattress is flat and free of any protruding foreign objects. Place the side with the product's "LOGO" facing upwards, with the power cord protruding from the side of the bed and plugged in.

Attention:

* This product is suitable for general mattresses, thick mattresses, nursing mattresses, foam mattresses, and palm fiber mattresses, etc., and its operating sensitivity is within a thickness range of 8-45cm.

* Do not place the monitoring pad on the sheets or lie directly on the pad.

Network pairing

1. Connect to power and wait for the Bluetooth blue light to stay on.
2. Turn on Bluetooth, Wi-Fi, and location services on your phone.
3. Open the Anvenus Plus app and press Connect.
4. Scan the QR code in the upper right corner of the monitoring pad to access the Wi-Fi configuration page.
5. Follow the on-screen instructions to select a Wi-Fi name, enter the password, and then press Connect.

* The network needs to be connected to a 2.4GHz Wi-Fi network; the configuration does not support the 5GHz band.

* Wi-Fi requires web authentication.

* If the network does not have a password, click Connect after selecting the network name.

6. After successful network pairing, the Bluetooth indicator light will turn off, and the Connection indicator light will flash green.

* Do not operate the device or sit/lie in bed during network pairing.


* Please complete network pairing within 1 minute, otherwise Bluetooth will automatically turn off. After Bluetooth turns off, you need to reconnect the power and wait for the Bluetooth blue light to stay on before re-pairing.

* If network pairing fails, please try reconnecting the power and pairing again.

* If network pairing fails, please check if Bluetooth and geolocation are enabled on your phone.

Place monitoring mattresses


1. Please contact your dealer for the login URL and account password.

2 After logging in, you will be taken to the main page, where you can view the monitoring screen to check the device usage status and user vital signs information.

3 Select the equipment and bed from the drop-down menu under the user name field in the management system.

Add a device and a bed, then bind it. Press "X" to return to the main page.

4. Click the bed card on the main page to enter the calibration interface. Calibration must be performed while out of bed. The PD value while out of bed is in the range of 850~950. During the calibration process, please do not operate or sit/lie on the bed. It is normal for vital signs data to appear during the calibration.

Precautions

1. If this product is not used for a long time, please disconnect the power supply.
2. This product is powered by a 5V-2A power plug (adapter). During installation, ensure proper grounding and avoid connecting it to the same power plug as high-power electrical appliances.
3. Because the sensor uses high-precision special optical fiber as a medium, it must be placed on a flat and stable surface for use. Please place bedding (mattresses) on top of this product. Placing the monitoring pad directly on hard objects or items with sharp protrusions may damage the core components of the product, resulting in inaccurate monitoring results.
4. When installing or moving the monitoring pad, Do not bent, rubbed, scratched, or torn, otherwise the product's lifespan will be affected.
5. The sensing part of the monitoring pad is designed without electricity, so it will not overheat and has no risk of fire, but it should still not be used in conjunction with flammable fire sources or high heat sources (such as electric blankets).
6. The surface of the monitoring pad is made of waterproof material, but it should be used in a dry environment and cannot be washed with water; if cleaning is necessary, the surface can be cleaned with a damp cloth.
7. The monitoring pad has high detection sensitivity. During use, please keep it away from mechanical vibration sources or noise sources (such as high-power speakers), otherwise it may damage the core components.
8. It should not be disposed of as household waste. The discarded equipment must be taken to an appropriate recycling center for disposal.
9. The monitoring pad is suitable for general mattresses, thick mattresses, nursing mattresses, sponge mattresses, and coir mattresses. The operating sensor is within the thickness range of 8-45cm. Do not place the pad on the bed sheet, and do not lie directly on the monitoring pad.
10. During monitoring, try to keep the chest area within the monitoring pad area and lie quietly flat.
11. If there are significant changes to the bedding or if the device is not accurately positioned for use, please recalibrate it.
12. Do not attempt to repair this product yourself. If the product malfunctions, please contact your dealer immediately.

Operation And Maintenance

1. Keep the product away from moisture. Excessive moisture or condensed Water vapor will cause damage to the parts of the product.
2. When in use, make sure there is no other source of vibration in bed.
3. Turn off or turn down the brightness of the AI monitor screen to enhance its life expectancy.
4. Keep the monitor in a well-ventilated area, away from excessive light, heat or moisture.
5. Clean the monitor screen gently with a clean cloth dampened with water. Do not spray any liquid directly onto the monitor screen.

Warranty Card

WARRANTY SERVICE CARD		
Customer name		Contact information
Product model		Warranty time
Purchase date		
Contact address		
Fault description		

*Please keep this card properly, and present this card if you need customer service.